

# Condensyn User Guide

## (Waterless Condensers)



- Utilising these waterless condensers can save the KCL over 1.5 million litres of water per year!
- Effective condensing is achieved through increased surface area combined with thicker glass

### Warnings and Guidance

- 2 units may be stacked on top of each other for longer term critical low loss experiments
- They are fragile – please use carefully and report any breakage to your local Safety Officer
- They come with two neck sizes – B19 and B24. Seek appropriate adapters if necessary
- Take caution after usage to not burn yourself as they can become hot
- It is recommended that to prevent breakthrough when using solvents of below 60 C boiling point that the heating bath/block temperature differential is kept below 10 C, e.g. for dichloromethane (bp 40 C) the heating bath/block temperature should not exceed 50 C. For diethyl ether this differential should be no more than 4 C

200mm CondenSyn is for reactions of 50ml and below

350mm CondenSyn is for reactions of 500ml and below

420mm CondenSyn is for reactions of 1000ml and below

- Note that a reactions that uses a combination of a solvent with a low boiling point and a low volume / long heating time, may be more suited to water-cooled condensers. Seek guidance if unsure