

多環芳烴標準品(PAHs)



特點與優勢:

經濟利益

- 即用型
- 廣泛的產品區間
- 可用於 GC, GC-MS, HPLC and LC-MS
- 符合 EPA 或 EU 分析方
- 法可制定

技術利益

- 按照 EPA 方法生產
- 產品一致性 - 可溯源，帶認證
- 可用於 EPA 500, 600 and 8000 系列方法
- COA 和 MSDS 可線上下載

這些產品基於重量/體積的基本原理進行稱量製備而成。溶質和 Reagecon 技術人員使用經過校準的器具制備。Reagecon 持有實驗室移液器的 ISO17025 認證

(INAB Ref: 265C)。由此產生的校準重量證書是按照 ISO/IEC 17025 的要求發佈的。然後，再利用重量法和密度法得出認證結果。Reagecon 公司的密度測量法是經過 ISO/IEC17025 認證的。每種標準液的濃度均利用氣相色譜-質譜聯用儀 (GC-MS) 進行驗證。並且，GC-MS 的校正是用經過 ISO Guide 34 認證的並與這些產品 PAH 濃度值類似的高純度 PAH 標準液完成的。通過與 NIST 的質譜庫進行比對來確定每個分析物的質譜。



Shannon Free Zone, Shannon, Co. Clare, Ireland
Tel: +353 81 472622 Fax: +353 81 472642
Email: sales@reagecon.ie
www.reagecon.com

CERTIFICATE OF GRIVAMETRIC PREPARATION

PRODUCT: Indeno (1,2,3-cd) pyrene @ 2,000 µg/ml in toluene
PRODUCT No.: REPAH115
MATRIX: Toluene
LOT NO.: PAH11516B1
DATE OF PREPARATION: 08th February 2016
EXPIRY DATE: 28th February 2019

PREPARATION OF STANDARD:

All analytical measuring devices and instrumentation have been pre-calibrated. The actual concentrations reported below are based on this preparation methodology and compound impurities.

Compound	Purity %	Labelled µg/mL	Actual µg/mL
Indeno (1,2,3-cd) pyrene	99.6	2000	1998 ±2.5 %

The reported uncertainty (±2.5%) due to weighing, volumetric preparation and homogeneity is calculated in compliance with EURACHEM/CEC/AC/ Guide: Quantifying Uncertainty in Analytical Measurements as ±2.5%.

Page 1 of 2

多環芳烴 (PAHs)混標

產品編號	描述	濃度與溶劑	包裝規格
REPAH001 (16 種混標)	benzo[k]fluoranthene 苯並(k)熒蒽 207-08-9	苯: 二氯甲烷	1ml
	acenaphthene 芴 83-32-9	2000mg/l	
	acenaphthylene 芴烯 208-96-8	EPA610 方法	
	anthracene 蒽 120-12-7		
	fluorene 芴 86-73-7		
	naphthalene 萘 91-20-3		
	phenanthrene 菲 85-01-8		
	benzo[a]anthracene 苯並[a]蒽 56-55-3		
	benzo[a]pyrene 苯並[a]芘 50-32-8		
	chrysene 屈 218-01-9		
	fluoranthene 熒蒽 206-44-0		
	indeno[1,2,3-cd]pyrene 茛並 (1,2,3-cd) 芘 193-39-5		
	pyrene 芘 129-00-0		
	benzo[b]fluoranthene 苯並(b)熒蒽 205-99-2		
	benzo[ghi]perylene 苯並[g,h,i]芘 191-24-2		
	dibenz[a,h]anthracene 二苯蒽 53-70-3		
產品編號	描述	濃度與溶劑	包裝規格
REPAH093 (16 種同位素 標記 PAH 混 標)	Acenaphthene-D10 芴-D10	100mg/l ,	1ml
	Acenaphthylene-D8 芴烯-D8	溶於正己烷	
	Anthracene-D10 蒽-D10	EPA610 方法	
	Benzo(a)anthracene-D12 苯並 (a) 蒽-D12		
	Benzo(b)fluoranthene-D12 苯並 (b) 熒蒽-D12		
	Benzo(k)fluoranthene-D12 苯並 (k) 熒蒽-D12		
	Benzo(g,h,i)perylene-D12 苯並 (g,h,i) 芘-D12		
	Benzo(a)pyrene-D12 苯並 (a) 芘-D12		
	Chrysene-D12 屈-D12		
	Dibenz(a,h)anthracene-D14 二苯並 (a,h) 蒽-D14		
	Fluoranthene-D10 熒蒽-D10		
	Fluorene-D10 芴-D10		
	Indeno(1,2,3-c,d)pyrene-D12 茛並[1,2,3-c,d]芘-D12		
	Naphthalene-D8 萘-D8		
	Phenanthrene-D10 菲-D10		
	Pyrene-D10 芘-D10		

多環芳烴(PAHs)混標

產品編號	描述	濃度 μg/ml	溶劑	包裝規格
REAPH003 (16 種混標)	benzo[k]fluoranthene 苯並(k)熒蒽 207-08-9	100	甲醇：丙酮	1ml
	acenaphthene 芴 83-32-9	1000		
	acenaphthylene 芴烯 208-96-8	2000		
	anthracene 蒽 120-12-7	100		
	fluorene 芴 86-73-7	200		
	naphthalene 萘 91-20-3	1000		
	phenanthrene 菲 85-01-8	100		
	benzo[a]anthracene 苯並[a]蒽 56-55-3	100		
	benzo[a]pyrene 苯並[a]芘 50-32-8	100		
	chrysene 屈 218-01-9	100		
	fluoranthene 熒蒽 206-44-0	200		
	indeno[1,2,3-cd]pyrene 茛並(1,2,3-cd)芘 193-39-5	100		
	pyrene 芘 129-00-0	100		
	benzo[b]fluoranthene 苯並(b)熒蒽 205-99-2	200		
	benzo[ghi]perylene 苯並[g,h,i]芘 191-24-2	200		
	dibenz[a,h]anthracene 二苯蒽 53-70-3	200		
產品編號	描述	濃度 μg/ml	溶劑	包裝規格
REPAH009 (18 種混標)	Acenaphthene 芴 83-32-9	1000	甲苯	1ml
	Anthracene 蒽 120-12-7			
	Benzo(a)anthracene 苯並[a]蒽 56-55-3			
	Chrysene 屈 218-01-9			
	Fluoroanthene 熒蒽 206-44-0			
	Fluorene 芴 86-73-7			
	Naphthalene 萘 91-20-3			
	Phenanthrene 菲 85-01-8			
	Pyrene 芘 129-00-0			
	Benzo(a)pyrene 苯並[a]芘 50-32-8			
	Benzo(b)fluoroanthene 苯並(b)熒蒽 205-99-2			
	Benzo(g,h,i)perylene 苯並[g,h,i]芘 191-24-2			
	Dibenzo(a,h)anthracene 二苯蒽 53-70-3			
	Benzo(k)fluoroanthene 苯並(k)熒蒽 207-08-9			
	Indeno(1,2,3-cd)pyrene 茛並(1,2,3-cd)芘 193-39-5			
	Acenaphthylene 芴烯 208-96-8			
	Benzo(j)fluoroanthene 苯並(j)熒蒽 205-82-3			
	Benzo(e)pyrene 苯並[e]芘 192-97-2			

多環芳烴(PAHs)單組分標準品

產品編號	描述	濃度 µg/ml	溶劑	包裝規格
REPAH3601	暈苯 191-07-1	100	甲苯	1ml
REPAH6701	1-甲基暈苯 13119-86-3	200	甲苯	1ml
REPAH3603	二苯並[a,j]暈苯 190-72-7	200	甲苯	1ml
REPAH3605	蔡基[2,3-a]暈苯 190-74-9	30	1,2,4-三氯苯	1ml
REPAH3701	1-硝基暈苯 81316-84-9	100	甲苯	1ml
REPAH3901	9-氯-9H-芴 6630-65-5	50	異辛烷	1ml
REPAH3902	2-氯蒽 17135-78-3	50	異辛烷	1ml
REPAH3903	9-氯菲 947-72-8	50	異辛烷	1ml
REPAH3904	6-氯苯[a]芘 21248-01-1	50	異辛烷	1ml
REPAH3905	1-氯芘 34244-14-9	50	異辛烷	1ml
REPAH3906	3-氯熗蒽 25911-51-7	50	異辛烷	1ml
REPAH3801	苯並蔥酮 82-05-3	1000	異辛烷	1ml
REPAH3804	異蔥酮紫 128-64-3	200	甲苯	1ml
REPAH3805	蔥酮紫 128-64-3	200	甲苯	1ml

多環芳烴(PAHs)單標

產品編號	描述	濃度 μg/ml	溶劑	包裝規格
REPAH101	芴 83-32-9	2000	甲苯	1ml
REPAH102	蒽 120-12-7	2000	甲苯	1ml
REPAH103	苯並[a]蒽 56-55-3	2000	甲苯	1ml
REPAH104	屈 218-01-9	2000	甲苯	1ml
REPAH105	芘 206-44-0	2000	甲苯	1ml
REPAH106	芴 86-73-7	2000	甲苯	1ml
REPAH107	蔡 91-20-3	2000	甲苯	1ml
REPAH108	菲 85-01-8	2000	甲苯	1ml
REPAH109	芘 129-00-0	2000	甲苯	1ml
REPAH110	苯並(a)芘 50-32-8	2000	甲苯	1ml
REPAH112	苯並[g,h,i]芘 191-24-2	2000	甲苯	1ml
REPAH113	二苯並(a,h)蒽 53-70-3	2000	甲苯	1ml
REPAH114	苯並(a)芘 50-32-8	100	乙腈	1ml
REPAH115	蒽並(1,2,3-cd)芘 193-39-5	2000	甲苯	1ml
REPAH116	芴烯 208-96-8	2000	甲苯	1ml
REPAH118	苯並[a]蒽 56-55-3	10	乙腈	1ml
REPAH119	苯並(a)芘 50-32-8	100	甲苯	1ml
REPAH150	2-乙醯氨基芴 53-96-3	1000	甲醇	1ml
REPAH151	2-乙醯氨基芴 53-96-3	2000	甲醇	1ml
REPAH152	7,12-二甲基苯並[a]蒽 57-97-6	1000	二氯甲烷:苯(1:1)	1ml
REPAH1141	苯並(k)芘 207-08-9	2000	甲苯	1ml
REPAH153	7,12-二甲基苯並[a]蒽 57-97-6	2000	二氯甲烷:苯(1:1)	1ml
REPAH134	蒽 95-13-6	2000	甲苯	1ml
REPAH154	二苯並(a,j)丫啶 224-42-0	1000	二氯甲烷	1ml
REPAH155	二苯並(a,j)丫啶 224-42-0	2000	二氯甲烷	1ml
REPAH156	二苯並[a,e]芘 192-65-4	1000	二氯甲烷:苯(1:1)	1ml
REPAH157	二苯並[a,e]芘 192-65-4	2000	二氯甲烷:苯(1:1)	1ml
REPAH1010	芘-d10 1718-52-1	200	異辛烷	1ml
REPAH147	芘 206-44-0	2000	甲苯	1ml
REPAH158	芘 206-44-0	1000	二氯甲烷	1ml
REPAH159	芘 206-44-0	2000	二氯甲烷	1ml
REPAH9001-I	2-氟-6-甲基蔡 324-42-5	100	異辛烷	1ml
REPAH9002-I	5-氟芴烯 17521-01-6	100	甲苯	1ml
REPAH9003-I	4-氟聯苯甲烷 587-79-1	100	甲苯	1ml
REPAH9004-I	2-氟芴 343-43-1	100	甲苯	1ml
REPAH9005-I	2-氟代二甲苯 3794-15-8	100	甲苯	1ml
REPAH9006-I	4,4'-二氟二苯甲烷 457-68-1	100	甲苯	1ml
REPAH9007-I	2-氟菲 523-41-1	100	甲苯	1ml
REPAH9008-I	3-氟菲 440-40-4	100	甲苯	1ml
REPAH9009-I	4-氟菲 521-66-4	100	甲苯	1ml
REPAH9010-I	3-氟-6-甲基菲 84194-32-2	50	異辛烷	1ml
REPAH1203	3-甲基膽蒽 56-49-5	100	乙腈	1ml

產品編號	描述	濃度 µg/ml	溶劑	包裝規格
REPAH9011-I	3-氟茱萸 1691-66-3	100	甲苯	1ml
REPAH9012-I	1-氟茱萸 1691-65-2	100	甲苯	1ml
REPAH9013-I	1-氟茱萸 1051914-84-1	100	甲苯	1ml
REPAH9014-I	3-氟茱萸 36288-22-9	100	甲苯	1ml
REPAH9015-I	9-氟-5-甲基茱萸 64977-46-8	50	異辛烷	1ml
REPAH9016-I	9-氟苯並(k)茱萸 113600-15-0	100	甲苯	1ml
REPAH7001	5-氟-3-甲基苯並[b]噻吩 17514-53-5	100	異辛烷	1ml
REPAH1009	芴-d10 81103-79-9	200	異辛烷	
REPAH7002	5-氟-2,3-二甲基苯並噻吩 1314962-98-5	100	異辛烷	1ml
REPAH7107	芴-d10 15067-26-2	1000	甲醇	1ml
REPAH7003	2-氟二苯噻吩 177586-38-8	100	甲苯	1ml
REPAH7101	1-甲基萘-d10 38072-94-5	1000	異辛烷	1ml
REPAH7102	2-甲基萘-d10 7297-45-2	1000	異辛烷	1ml
REPAH7103	1,8-二甲基萘-d12 104489-29-4	1000	異辛烷	1ml
REPAH7104	2,6-二甲基萘-d12 350820-12-1	1000	異辛烷	1ml
REPAH7105	9-甲基萸-d12 6406-97-9	1000	異辛烷	1ml
REPAH7106	1-甲基茱萸-d9 210487-07-3	100	甲苯	1ml
REPAH1102	三亞苯-d12 17777-56-9	1000	甲苯	1ml
REPAH1103	苯並(e)茱萸-d12 205440-82-0	100	甲苯	1ml
REPAH1104	苯並(b)茱萸-d12 93951-98-5	100	甲苯	1ml
REPAH1105	苯並(k)茱萸-d12 93952-01-3	100	甲苯	1ml
REPAH1106	苯並(ghi)茱萸-d12 93951-66-7	100	甲苯	1ml
REPAH1107	苯並(ghi)茱萸-d12 93951-66-7	200	甲苯	1ml
REPAH1108	茛並[1,2,3-CD]茱萸-d12 203578-33-0	100	甲苯	1ml
REPAH1109	二苯並(a,h)萸-d14 13250-98-1	100	異辛烷	1ml
REPAH1110	二苯並(a,h)萸-d14 13250-98-1	100	甲苯	1ml
REPAH1111	暈苯-d12 16083-32-2	100	甲苯	1ml
REPAH1112	二苯並(a,i)茱萸-d14 158776-07-9	100	甲苯	1ml
REPAH1113	聯苯-d10 1486-01-7	1000	甲苯	1ml
REPAH1114	鄰三聯苯-d14 5142-76-6	100	甲苯	1ml
REPAH1115	間三聯苯-d14 17714-84-0	100	甲苯	1ml
REPAH1116	對三聯苯-d14 1718-51-0	100	甲苯	1ml
REPAH1117	對三聯苯-d14 1718-51-0	1000	甲苯	1ml
REPAH1118	2,2'-聯萘-d14 210487-05-1	100	甲苯	1ml
REPAH1201	噻唑-d8 97960-57-1	1000	甲苯	1ml
REPAH1202	吡啶-d9 34749-75-2	1000	甲苯	1ml
REPAH1301	1-硝基萘-d7 80789-77-1	1000	甲苯	1ml
REPAH1302	2-甲基-1-硝基萘-d9	100	異辛烷	1ml
REPAH1303	2-硝基芴-d9 128008-87-7	100	甲苯	1ml
REPAH1304	2-硝基芴-d9 128008-87-7	1000	甲苯	1ml

多環芳烴(PAHs)單組分標準品

產品編號	描述	濃度 µg/ml	溶劑	包裝規格
REPAH1305	9-硝基菲-d9	10	甲苯	1ml
REPAH1306	9-硝基蒽-d9 220381-38-4	100	甲苯	1ml
REPAH1307	1-硝基芘-d9 93487-20-8	100	甲苯	1ml
REPAH1308	3-硝基熒蒽-d9 350820-11-0	100	甲苯	1ml
REPAH1309	1-硝基三亞苯-d11 350820-11-0	100	異辛烷	1ml
REPAH1310	6-硝基蒽-d11 203805-92-9	100	甲苯	1ml
REPAH1311	6-硝基苯並(a)芘-d11 352431-12-0	100	甲苯	1ml
REPAH1401	1-氨基萸-d7 78832-53-8	1000	甲苯	1ml
REPAH1402	2-氨基萸-d7 93951-94-1	1000	異辛烷	1ml
REPAH1001	1-硝基萸 86-57-7	100	甲苯	1ml
REPAH1002	2-硝基萸 581-89-5	100	甲苯	1ml
REPAH1003	1-甲基-4-硝基萸 880-93-9	100	甲醇	1ml
REPAH1004	1-甲基-5-硝基萸 91137-27-8	100	甲醇	1ml
REPAH1005	1-甲基-6-硝基萸 105752-67-8	100	甲醇	1ml
REPAH1006	2-甲基-1-硝基萸 881-03-8	100	甲醇	1ml
REPAH1007	2-甲基-4-硝基萸 13615-38-8	100	甲醇	1ml
REPAH1008	1,5-二硝基萸 605-71-0	100	甲苯	1ml
REPAH5001	1-甲基熒蒽 1730-37-6	1000	甲苯	1ml
REPAH5002	2-甲基熒蒽 1430-97-3	1000	甲苯	1ml
REPAH5003	4-甲基熒蒽 1556-99-6	1000	異辛烷	1ml
REPAH5004	9-甲基熒蒽 2523-37-7	1000	異辛烷	1ml
REPAH5005	1,7-二甲基熒蒽 442-66-0	500	異辛烷	1ml
REPAH5006	9-乙基熒蒽 2294-82-8	1000	異辛烷	1ml
REPAH5007	9-n-丙基熒蒽 4037-45-0	1000	異辛烷	1ml
REPAH5008	9-n-丁基熒蒽 3952-42-9	1000	異辛烷	1ml
REPAH5009	9,9-二辛基熒蒽 123863-99-0	1000	異辛烷	1ml
REPAH5010	9,9'-二亞甲基芴 746-47-4	1000	甲苯	1ml
REPAH4101	11H-苯並(a)芴 238-84-6	1000	甲苯	1ml
REPAH4103	11H-苯並(b)芴 243-17-4	200	甲苯	1ml
REPAH4104	7H-苯並(c)芴 205-12-9	200	甲苯	1ml
REPAH4105	9-苯基芴 789-24-2	1000	異辛烷	1ml
REPAH4201	2-硝基芴 607-57-8	100	甲苯	1ml
REPAH4202	2,7-二硝基芴 5405-53-8	100	甲苯	1ml
REPAH4203	2-硝基-9-芴酮 3096-52-4	100	甲苯	1ml
REPAH4401	菲 85-01-8	1000	異辛烷	1ml
REPAH5101	1-甲基菲 832-69-9	1000	異辛烷	1ml
REPAH5102	2-甲基菲 2531-84-2	1000	異辛烷	1ml
REPAH5103	3-甲基菲 832-71-3	1000	異辛烷	1ml
REPAH5104	4-甲基菲 832-64-4	500	異辛烷	1ml
REPAH5105	9-甲基菲 883-20-5	1000	異辛烷	1ml

產品編號	描述	濃度 µg/ml	溶劑	包裝規格
REPAH5201	1,2-二甲基菲 20291-72-9	500	異辛烷	1ml
REPAH5202	1,3-二甲基菲 16664-45-2	500	異辛烷	1ml
REPAH5203	1,4-二甲基菲 22349-59-3	500	異辛烷	1ml
REPAH5204	1,5-二甲基菲 66271-87-2	50	異辛烷	1ml
REPAH5205	1,6-二甲基菲 20291-74-1	500	異辛烷	1ml
REPAH5206	1,7-二甲基菲 483-87-4	500	異辛烷	1ml
REPAH5207	1,8-二甲基菲 7372-87-4	500	異辛烷	1ml
REPAH5208	1,9-二甲基菲 20291-73-0	500	異辛烷	1ml
REPAH5209	2,3-二甲基菲 3674-65-5	50	異辛烷	1ml
REPAH5210	2,4-二甲基菲 15254-64-5	500	異辛烷	1ml
REPAH5211	2,5-二甲基菲 3674-66-6	50	異辛烷	1ml
REPAH5212	2,6-二甲基菲 17980-16-4	50	異辛烷	1ml
REPAH5213	2,7-二甲基菲 1576-69-8	500	異辛烷	1ml
REPAH5214	2,9-二甲基菲 17980-09-5	50	異辛烷	1ml
REPAH5215	2,10-二甲基菲 2479-54-3	50	異辛烷	1ml
REPAH5216	3,4-二甲基菲 66291-31-4	50	異辛烷	1ml
REPAH5217	3,5-二甲基菲 33954-06-2	50	異辛烷	1ml
REPAH5218	3,6-二甲基菲 1576-67-6	500	異辛烷	1ml
REPAH5219	3,9-二甲基菲 66291-32-5	500	異辛烷	1ml
REPAH5220	3,10-二甲基菲 66291-33-6	500	異辛烷	1ml
REPAH5221	4,9-二甲基菲 66291-34-7	50	異辛烷	1ml
REPAH5222	4,10-二甲基菲 23189-63-1	50	異辛烷	1ml
REPAH5223	9,10-二甲基菲 604-83-1	500	異辛烷	1ml
REPAH5224	3-乙基菲 1576-68-7	500	異辛烷	1ml
REPAH5225	9-乙基菲 3674-75-7	1000	異辛烷	1ml
REPAH5301	1,2,4-三甲基菲 23189-64-2	200	異辛烷	1ml
REPAH5302	1,2,5-三甲基菲 146448-87-5	50	異辛烷	1ml
REPAH5303	1,2,7-三甲基菲 60490-68-8	50	異辛烷	1ml
REPAH5304	1,2,6-三甲基菲 30436-55-6	500	異辛烷	1ml
REPAH5305	1,2,8-三甲基菲 20291-75-2	500	異辛烷	1ml
REPAH5306	1,2,9-三甲基菲 146448-88-6	500	異辛烷	1ml
REPAH5307	1,3,4-三甲基菲 66271-45-2	200	異辛烷	1ml
REPAH5308	2,6,9-三甲基菲 66271-32-7	200	異辛烷	1ml
REPAH5309	2,6,9-三甲基菲 66271-32-7	1000	異辛烷	1ml
REPAH5310	9-n-丙基菲 17024-03-2	1000	異辛烷	1ml
REPAH5401	1,2,6,9-四甲基菲 204256-39-3	500	異辛烷	1ml
REPAH5402	9-n-丁基菲 10394-57-7	500	異辛烷	1ml
REPAH5403	蒽烯 483-65-8	500	異辛烷	1ml
REPAH5404	1,9-甲基-5-乙基菲	50	異辛烷	1ml

多環芳烴(PAHs)單組分標準品

產品編號	描述	濃度 µg/ml	溶劑	包裝規格
REPAH5405	1,9-二甲基-7-乙基菲 46870-52-4	50	異辛烷	1ml
REPAH1501	4H-環戊二烯[def]菲 203-64-5	500	異辛烷	1ml
REPAH1502	1H-環戊二烯[i]菲 235-92-7	500	異辛烷	1ml
REPAH1503	苯並[c]菲 195-19-7	200	甲苯	1ml
REPAH1504	2-甲基環戊二烯[i]菲 121254-39-5	500	異辛烷	1ml
REPAH1505	三亞苯 217-59-4	200	甲苯	1ml
REPAH1506	3-甲基菲並[3,4-C]菲 83844-21-7	100	甲苯	1ml
REPAH1701	1-甲氧基菲 834-99-1	1000	異辛烷	1ml
REPAH1702	2-甲氧基菲 13837-48-4	1000	異辛烷	1ml
REPAH1703	3-甲氧基菲 835-06-3	1000	異辛烷	1ml
REPAH1704	4-甲氧基菲 15638-06-9	1000	異辛烷	1ml
REPAH1705	9-甲氧基菲 5085-74-5	1000	異辛烷	1ml
REPAH1801	2-硝基菲 17024-18-9	1000	異辛烷	1ml
REPAH1802	3-硝基菲 17024-19-0	300	異辛烷	1ml
REPAH1803	9-硝基菲 954-46-1	1000	異辛烷	1ml
REPAH1804	5-硝基苯[c]菲 64356-30-5	100	甲苯	1ml
REPAH1805	1-硝基三亞苯	50	異辛烷	1ml
REPAH1901	蒽 120-12-7	1000	異辛烷	1ml
REPAH5501	1-甲基蒽 610-48-0	1000	異辛烷	1ml
REPAH5502	2-甲基蒽 613-12-7	1000	異辛烷	1ml
REPAH5504	9-甲基蒽 779-02-2	1000	異辛烷	1ml
REPAH5601	1,2-二甲基蒽 53666-94-7	200	甲苯	1ml
REPAH5602	1,3-二甲基蒽 610-46-8	200	甲苯	1ml
REPAH5603	1,4-二甲基蒽 781-92-0	200	甲苯	1ml
REPAH5604	1,5-二甲基蒽 15815-48-2	200	甲苯	1ml
REPAH5605	2,3-二甲基蒽 613-06-9	200	甲苯	1ml
REPAH5606	2,7-二甲基蒽 782-23-0	200	甲苯	1ml
REPAH5607	9,10-二甲基蒽 781-43-1	200	甲苯	1ml
REPAH5608	2-乙基蒽 52251-71-5	1000	異辛烷	1ml
REPAH5701	1,2,4-三甲基蒽 20153-28-0	200	甲苯	1ml
REPAH5702	1,2,3,4-四甲基蒽 66553-01-3	200	甲苯	1ml
REPAH5703	2,3,6,7-四甲基蒽 15254-25-8	200	甲苯	1ml
REPAH5704	2,3,9,10-四甲基蒽 66552-77-0	200	甲苯	1ml
REPAH5705	2-叔丁基蒽 18801-00-8	1000	異辛烷	1ml
REPAH5801	1-甲基苯並[a]蒽 2498-77-3	50	甲苯	1ml
REPAH5802	5-甲基苯並[a]蒽 2319-96-2	200	甲苯	1ml
REPAH5803	6-甲基苯並[a]蒽 316-14-3	200	甲苯	1ml

產品編號	描述	濃度 µg/ml	溶劑	包裝規格
REPAH5804	7-甲苯並[a]蒽 2541-69-7	50	甲苯	1ml
REPAH5805	10-甲苯並[a]蒽 2381-15-9	200	甲苯	1ml
REPAH5806	3,9-二甲苯並[a]蒽 316-51-8	200	甲苯	1ml
REPAH5807	7,12-二甲苯並[a]蒽 57-97-6	200	甲苯	1ml
REPAH2001	二苯並[a,c]蒽 215-58-7	100	甲苯	1ml
REPAH2002	四苯並[a,c,h,i]蒽 215-11-2	200	甲苯	1ml
REPAH2101	2-硝基蒽 3586-69-4	200	甲苯	1ml
REPAH2102	9-硝基蒽 602-60-8	100	甲苯	1ml
REPAH2103	9-甲基-10-硝基蒽 84457-22-7	100	甲苯	1ml
REPAH2104	9,10-二硝基蒽 33685-60-8	100	甲苯	1ml
REPAH2201	7-硝基苯並[a]蒽 20268-51-3	100	甲苯	1ml
REPAH2202	7-硝基二苯並[a,h]蒽 63041-91-8	100	甲苯	1ml
REPAH2301	1,2,3,10b-Tetrahydrofluoranthene 20279-21-4	200	甲苯	1ml
REPAH5901	1-甲基熹蒽 25889-60-5	200	甲苯	1ml
REPAH5902	2-甲基熹蒽 33543-31-6	200	甲苯	1ml
REPAH5903	3-甲基熹蒽 1706-01-0	200	甲苯	1ml
REPAH5904	3-乙基熹蒽 20496-16-6	200	甲苯	1ml
REPAH2401	苯並[a]熹蒽 203-33-8	200	甲苯	1ml
REPAH2402	苯並[ghi]熹蒽 203-12-3	200	甲苯	1ml
REPAH2403	2-苯基熹蒽 4914-78-7	200	甲苯	1ml
REPAH2502	二苯並[a,e]熹蒽 5385-75-1	200	甲苯	1ml
REPAH2506	蒽並[1,2,3-cd]熹蒽 193-43-1	200	甲苯	1ml
REPAH2601	蔡並[1,2-b]熹蒽 111189-32-3	200	甲苯	1ml
REPAH2602	蔡並[1,2-k]熹蒽 238-04-0	200	甲苯	1ml
REPAH2603	蔡並[2,3-b]熹蒽 206-06-4	200	甲苯	1ml
REPAH2604	蔡並[2,3-j]熹蒽 205-83-4	200	甲苯	1ml
REPAH2605	蔡並[2,3-k]熹蒽 207-18-1	200	甲苯	1ml
REPAH2701	1-硝基熹蒽 13177-28-1	100	甲苯	1ml
REPAH2702	2-硝基熹蒽 13177-29-2	100	甲苯	1ml
REPAH2703	3-硝基熹蒽 892-21-7	100	甲苯	1ml
REPAH6001	1-甲基芘 2381-21-7	200	甲苯	1ml
REPAH6002	4-甲基芘 3353-12-6	200	甲苯	1ml
REPAH6003	4,5-二甲基芘 15679-25-1	200	甲苯	1ml
REPAH6004	2,7-二甲基芘 15679-24-0	200	甲苯	1ml
REPAH6005	1-乙基芘 17088-22-1	200	甲苯	1ml
REPAH6006	1-n-丙基芘 42211-33-6	1000	甲苯	1ml
REPAH6007	1-n-丁基芘 35980-18-8	200	甲苯	1ml
REPAH6101	6-甲基苯[a]芘 2381-39-7	200	甲苯	1ml
REPAH6102	7-甲基苯[a]芘 63041-77-0	200	甲苯	1ml
REPAH6103	7,10-二甲基苯[a]芘 63104-33-6	50	甲苯	1ml
REPAH2801	二苯並[a,e]芘 192-65-4	200	甲苯	1ml
REPAH2802	二苯並[a,h]芘 189-64-0	200	甲苯	1ml
REPAH2803	二苯並[a,i]芘 189-55-9	200	甲苯	1ml

產品編號	描述	濃度 µg/ml	溶劑	包裝規格
REPAH2804	二苯並[a,l]芘 191-30-0	200	甲苯	1ml
REPAH2805	二苯並[e,l]芘 192-51-8	200	甲苯	1ml
REPAH2901	環戊烯(c,d)芘 27208-37-3	50	甲苯	1ml
REPAH2902	蔡[2,3-a]芘 196-42-9	200	甲苯	1ml
REPAH2903	蔡[2,3-e]芘 193-09-9	200	甲苯	1ml
REPAH2904	2,3-Peri-亞蔡基-芘 188-89-6	200	甲苯	1ml
REPAH2905	2,3,7,8-二-(peri-亞蔡基)-芘 362052-03-7	200	甲苯	1ml
REPAH3001	1-羥基芘 5315-79-7	200	甲苯	1ml
REPAH3003	3-羥基苯[a]芘 13345-21-6	50	甲苯	1ml
REPAH3101	1-硝基芘 5522-43-0	100	甲苯	1ml
REPAH3102	2-硝基芘 789-07-1	100	甲苯	1ml
REPAH3103	4-硝基芘 57835-92-4	100	甲苯	1ml
REPAH3104	1,3-二硝基芘 75321-20-9	100	甲苯	1ml
REPAH3105	1,6-二硝基芘 42397-64-8	100	甲苯	1ml
REPAH3106	1,8-二硝基芘 42397-65-9	100	甲苯	1ml
REPAH6201	1-甲基蒽 3351-28-8	200	甲苯	1ml
REPAH6202	2-甲基蒽 3351-32-4	200	甲苯	1ml
REPAH6203	3-甲基蒽 3351-31-3	200	甲苯	1ml
REPAH6204	4-甲基蒽 3351-30-2	200	甲苯	1ml
REPAH6205	5-甲基蒽 3697-24-3	200	甲苯	1ml
REPAH6206	6-甲基蒽 1705-85-7	200	甲苯	1ml
REPAH6301	6-乙基蒽 2732-58-3	1000	甲苯	1ml
REPAH6302	1,3,6-三甲基蒽	1000	異辛烷	1ml
REPAH6303	6-n-丙基蒽 6910-41-4	1000	甲苯	1ml
REPAH6304	6-n-丁基蒽 6901-71-9	1000	甲苯	1ml
REPAH3201	苯並[a]蒽 213-46-7	100	甲苯	1ml
REPAH3202	苯並[b]蒽 214-17-5	200	甲苯	1ml
REPAH3203	苯並[c]蒽 194-69-4	200	甲苯	1ml
REPAH3204	苯並[g]蒽 196-78-1	200	甲苯	1ml
REPAH3205	二苯並[g,p]蒽 191-68-4	200	甲苯	1ml
REPAH3301	Anthanthrene 蒽嵌蒽 191-26-4	200	甲苯	1ml
REPAH6401	6-甲基蒽嵌蒽 31927-64-7	200	甲苯	1ml
REPAH3302	6-硝基蒽嵌蒽 111438-94-9	100	甲苯	1ml
REPAH6501	1-n-己基芘 143076-98-6	200	甲苯	1ml
REPAH3401	苯並[b]芘 197-70-6	200	甲苯	1ml
REPAH3402	二苯並[b,ghi]芘 5869-30-7	200	甲苯	1ml
REPAH3403	二苯並[e,ghi]芘 190-95-4	200	甲苯	1ml
REPAH3404	蔡並[8,1,2-bcd]芘 88-89-6	200	甲苯	1ml
REPAH3405	蔡並[1,2,3,4-ghi]芘 190-84-1	200	甲苯	1ml
REPAH3501	1-硝基芘 35337-20-3	100	甲苯	1ml
REPAH3502	3-硝基芘 20589-63-3	100	甲苯	1ml
REPAH6601	2,9-二甲基芘 1679-02-3	100	甲苯	1ml